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- The Bebeconfort Mini multifunction thermometer is to be operated by consumers in the home setting and primary care setting as screening tool. This manual, accessories, Directions for Use, all precautionary information, and specifications should be read before use.
- The thermometer is intended to measure the temperature of human body, at the surface of the forehead skin or ear, of people with any age and any weight (except of any other medical advice). The patient can use the thermometer by himself.
The measurement mode is direct.
Do not use it for any other purpose.
- This product is intended in the home setting and primary care setting as screening tool.
- Do not use the thermometer if it malfunctions or has been damaged in any matter.
- When the ambient temperature of the thermometer changes too much, such as moving the Thermometer from one place of lower temperature to another place of higher temperature, allow the thermometer to remain in a room for 30 minutes where the temperature is between 10°C to 40°C (50°F- 104°F).
- Remove primary batteries if equipment is not likely to be used for long time.
- This product is not waterproof, do not be immersed in water or other liquid; If cleaning and disinfection, please follow the «Cleaning, care and storage» section requirements.
- Do not touch the sensor of infrared detection with your fingers.
- If a cold compress on the forehead fever patients, or take other measures to cool down the temperature data will low, should be avoided in this case to measure body temperature.
- This product must be operated in a stable environment. If ambient environment changes significantly, check if there is fog on the sensor. If any, remove it in accordance with «Cleaning, Care and Storage» section.
- Do not near strong electrostatic field or strong magnetic fields, in order to avoid to affect the accuracy of the measurement data.
- The measuring result of this product is only for your reference. If you have any doubt, please measure the temperature in other methods.
- It may affect the accuracy of measurements when the forehead is covered by hair, perspiration, cap or scarf. Also the wax of ear and other thin layers will low the temperature. So please be sure use correct way in case of accuracy of measurements.
- Some people are not able to use ear mode properly for oily ears or for other reasons.
- If measure human forehead temperature , please select «body» mode; for measure other objects, liquids, food and other temperature please select «surface» mode.

WARNINGS

- The device should be kept out of the reach of children/pets. When not in use, store the device in a dry room and protect against extreme moisture, heat, lint, dust and direct sunlight. Never place any heavy objects on the device.
- Do not throw batteries into fire.
- Only use recommended batteries. Do not use rechargeable batteries.
- This thermometer does not replace the diagnostic of a medical professional.
- Do not fall, disassemble or modify the device.
- Do not use this device if you think it is damaged or notice anything unusual.
- This device comprises sensitive components and must be treated with caution. Observe the storage and operating conditions described in the 'Technical Specifications' section.
- Not servicing/maintenance while the Thermometer is in use.
- When using, shall not touch battery and the patient simultaneously.
- Do not use the device if it is damaged/degraded/loosened in any way.
- The continuous use of a damaged unit may cause injury, improper results, or serious danger.
- Based on the current science and technology, other potential allergic reactions are unknown.

- This equipment needs to be installed and put into service in accordance with the information provided in the user manual.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the thermometer, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
- If the use location is near (e.g. less than 1.5 km from) AM, FM or TV broadcast antennas, before using this equipment, it should be observed to verify that it is operating normally to assure that the equipment remains safe with regard to electromagnetic disturbances throughout the expected service life.

1- Overview

The Bebeconfort Mini multifunction thermometer has been designed to offer 3 operating modes of temperatures measurement: forehead, ear, surface and object. This thermometer uses the new generation of infrared scanning technology with contactless scanning, high measurement accuracy and ultra-speed measurement in 2 secs. This thermometer can be used to measure baby and adult body temperatures.

Description on Controls, Indicators, and Symbols in the home setting and as primary care setting screening tool


Figure 1: Structure description

1. Passive infrared probe
2. Cap
3. Battery cover
4. Liquid crystal display (LCD)
5. On/Memo button
6. Scan button



Figure 2: LCD description

1. Forehead mode
2. Ear mode
3. Memory indicator
4. Speaker
5. Battery indicator Data indicator
6. Measurement unit
7. Indicator of measurement result

Thermometer applications

Thermometer Model Number	Thermometer Style	Adult		Pediatric	
		Ear	Forehead	Ear	Forehead
HTD8219EU	Mini multifunction thermometer	√	√	√	√

SYMBOL EXPLANATION

	Batch Code		Manufacturer
	Consult Accompanying Documents		Type BF Applied Part
	Operating Temperature		Operating Humidity
	Storage Temperature		Storage Humidity
	Medical Device		Storage and Operating atmospheric pressure
	European Authorized Representative		Disposal of this product and used batteries should be carried out in accordance with the national regulations for the disposal of electronic products.
	Indicates this device is in compliance with Medical Device Directive 93/42/EEC. 0598 is the Notified Body Number		
	IP22 - The first number 2 : Protected against solid foreign objects of 12,5mm Ø.and greater. The second number 2 : Protected against vertically falling water drops when enclosure tilted up to 15°		
	Serial number Noted SN is printed inside the thermometer		

TECHNICAL SPECIFICATIONS	
Measurement Unit	°C
Operating mode	Adjusted mode (Forehead mode/ Ear mode) Direct mode (surface mode)
Measuring site	Forehead / Ear / Surface / Object
Reference Body Site	Oral
Rated output range	Body mode: 34.0°C~43.0°C Surface mode: 0°C~100.0°C
Laboratory Accuracy	Body mode: 34.0°C~34.9°C:±0.3°C; 35.0°C~42.0°C:±0.2°C; 42.1°C~43.0°C:±0.3°C; Surface mode: ±2°C
Display Resolution	0.1°C
Alarm of Background light	Green: 34-37.4°C Normal Temperature Yellow: 37.5°C-38.5°C Slight Fever Yellow with flash: 38.6-43.0°C High Fever Note: Surface mode is always with Green backlight.
Auto Power Off time	≤35s
Measuring time	Ear mode : ≤2S Forehead mode : ≤5S (for forehead scan mode) Surface mode : ≤5S
Memories	50

Power Supply Requirements	
Batteries	DC3.0V CR2032×1
Environmental	
Operating Condition	Operating Temperature:10°C~40°C, Relative Humidity≤85%, atmospheric pressure:70-106Kpa
Transport and Storage Condition	Storage Temperature:-20-55°C, Relative Humidity≤95%, atmospheric pressure:70-106Kpa

Accessory name	number
CR2032	1
User manual	1

Dimension and Weighting	
Weight (without batteries)	52g
Size	L:93.6mm X W:51.3mm X H:24.5mm

Compliance	
Item	Compliant with
Equipment classification	Safety Standards: IEC 60601-1:2005+A1:2012+A2:2020, IEC 60601-1-2:2014+A1:2020
Type of protection	Internally powered equipment (on battery power)
Front panel and case labeling	ISO 15223-1:2021
Temperature	ISO 80601-2-56:2017+A1:2018
Home healthcare environment	IEC 60601-1-11:2015+A1:2020

Calculated values of the indicators according to ISO 80601-2-56

Ear mode:

Indicators	Group A1	Group A2	Group B	Group C
Clinical bias, Δcb	0.110	0.105	0.100	0.110
Standard deviation, σ_j	0.133	0.110	0.120	0.133
Limits of agreement, LA	0.267	0.220	0.239	0.266
Clinical repeatability (for All Group), σ_f	0.121			

Forehead mode:

Indicators	Group A1	Group A2	Group B	Group C
Clinical bias, Δcb	0.170	0.160	0.150	0.154
Standard deviation, σ_j	0.108	0.099	0.086	0.146
Limits of agreement, LA	0.216	0.199	0.173	0.292
Clinical repeatability (for All Group), σ_f	0.108			

Note: the above value is calculated from clinical data of HTD8222EU.

Group A1: 0 months up to 3 months, Group A2: 3 months up to 1 year,

Group B: Older than 1 and younger than 5 years, Group C: Older than 5 years.

SAFETY CLASSIFICATION OF MEDICAL ELECTRICAL (ME) EQUIPMENT

Protection against electric shock	Internally powered ME equipment
Degree of protection	Type BF - Applied part(the probe)
Protection against harmful ingress of water or particulate matter	IP22
Mode of operation	Continuous operation

Note: Not intended to be sterilized. Not for use in an OXYGEN RICH ENVIRONMENT

2 - TEMPERATURE MESURING METHODS

1) Forehead scan mode

- 1- Install cap on the probe and turn on the thermometer by pressing the On/Memo button. The thermometer will perform self-test with all segments displayed for 1 second. See figure 1
- 2- After self-test is completed, the device is ready to measurement with " - - - ", see figure 2

Figure 1


Figure 2

Forehead temperature measurement :

Apply thermometer to the forehead and press the Scan button around 5 seconds from left to right or from right to left. A short "bip" will be heard when measurement is completed. Then the measured temperature will be displayed on the LCD screen with the appropriate backlight colour.

Backlight indications:

Green: 34-37.4°C Normal Temperature

Yellow: 37.5°C-38.5°C Slight Fever

Yellow with flash: 38.6-43.0°C High Fever

Note that both yellow and yellow with flash are warning.

2) Ear mode

- 1- Remove the cap from the probe and turn on the thermometer by pressing the On/Memo button. The thermometer will perform self-test with all segments displayed. See figure 3
- 2- After self-test is completed, the device is ready to measurement with " - - - ", see figure 4

Figure 3


Figure 4

Ear temperature measurement :

Insert the thermometer probe slowly into the ear canal and press the Scan button. A short "bip" will be heard when measurement is completed. Then the measured temperature will be displayed on the LCD screen with the appropriate backlight colour.

Backlight indications:

Green: 34-37.4°C Normal Temperature

Yellow: 37.5°C-38.5°C Slight Fever

Yellow with flash: 38.6-43.0°C High Fever

Note that both yellow and yellow with flash are warning.

Note: If the person for ear temperature measurement is more than one year old, please pull the ear up and back gently. If less than one year old, pull ear back gently. See figure 5


Figure 5

3) Surface/Object mode

1. Install cap on the probe and turn on the thermometer.

2. When the device is powered on, press and hold the On/Memo button and Scan button at the same time to cycle conversion between «body» mode and «surface» mode. For the display of surface mode, see Figure 6.

«body» mode is used for measuring human body temperature, the «surface» mode is used to measure the object's surface temperature. (The factory default is «body» mode).

Figure 6

Surface or object measurement:

Point the thermometer wherever you want to measure the temperature (in the middle of a room, in front of a plate of food, or a warm bath) and press the Scan button around 5s from left to right or from right to left. And a short «bip» will be heard, when measurement is completed. Then the measured temperature will be displayed on the LCD screen.

Note: Surface mode is always with Green backlight.

General notes:

- Please keep the forehead is free of hair, perspiration, cap, scarf and etc when forehead temperature measurement.
- Please keep the ear canal and probe surface clean when ear temperature measurement.
- People body should at least remain in the test environment more than 5 minutes, and measure after the ambient temperature is in accordance with the environment.
- If a cool compress or other measures were taken to cool down forehead of fever patient, the temperature data will low.
- Surround ambient temperature should be stable, keep away from where airflow is strong, such as fan, air conditioner vents etc.
- When the thermometer moving from one place of lower temperature to another place of higher temperature, should at least remain in the test environment more than 30 minutes, to be consistent with the ambient temperature .
- You cannot use the thermometer in place where the sun is strong.
- Recommended test at least 3 times and select the result that appears most.
- If measure body temperature, please select «ear» mode or «forehead» mode.
- If for some reason caused the forehead temperature is low, you can try to measure behind the ear or take ear mode to measure.

3 - RECALLING AND ERASING MEMORY DATA

The last temperature taken before the thermometer powers off is stored in memory.

- 1) In the boot state, short press the On/Memo button to view the history of measured values.
- 2) An empty memory cell shows "---°C".
- 3) Temperature readings can be stored in memory. Up to 50 temperature readings can be stored into the memory cells and automatically overwrite historical data.
- 4) In boot mode, press and hold the On/Memo and Scan button simultaneously until the LCD display «CLR», which means that all stored data is cleared completely.

4 - HOW TO SET SOUND MUTE OR NOT

In boot mode, press and hold On/Memo to set sound mute or not.

5 - TROUBLESHOOTING

MESSAGE	SITUATION	SOLUTION
	Temperature taken in is higher than Typical human temperature.	detected target temperature higher than 43.0°C
	Temperature taken in is lower than Typical human temperature.	detected target temperature lower than 34.0°C
	Operating temperature exceeds the range of specified temperature.	Environment temperature higher than 40°C or lower than 10°C
	Wrong function indicate	System failure found when System auto-check, the screen flicker, and automatic turn off.
	Battery capacity is too low. Taking Temperature is not allowed.	Install a new battery
	Error warning	The interval between the two measurements is not up to 5S
	Ambient temperature is unstable.	Ambient temperature is unstable. Please allow it to acclimate for at least 30 minutes before using it.
	Blank screen indicate	(1) Press ON button again. (2) Check the battery polarity. (3) Replace with a new battery.

6 - REPLACING THE BATTERY

- 1) Open and release battery cover following indicator using screw driver. Before changing the battery being sure the system is already power off.
- 2) Remove the battery and replace with a new one, CR2032, make sure align properly as indicated inside the battery compartment.
- 3) Slide the battery cover back in and lock the screw back using screw driver.
- 4) Do not dispose of used batteries in household waste. Take them to special local collection sites.
- 5) If system is locked after changing battery. You may not follow up the process of rule one. Just take off battery, waiting 30 sec, then load battery again.

WARNING!

Do not recharge, disassemble or dispose of in fire.

7 - CLEANING, CARE AND STORAGE

The lens is very delicate. It is very important to protect the lens from dirt and damage. The Infrared probe is very delicate. Do not touch or use tools press it. Must be carefully protected otherwise it will affect the accuracy of the measure.

Use a clean, soft cloth to clean the surface of the device and LCD. Do not use solvents or immerse the device into water or other liquids.

It is recommended to store the thermometer in a dry location free from dust.

Always keep the thermometer within the storage temperature range (- 20°C to 55°C) and humidity range ($\leq 95\%$ non-condensing). At least 30 min required of equipment to warm from the minimum storage temperature between uses until it is ready for intended use.

At least 30 min required for me equipment to cool from the maximum storage temperature between uses until it is ready for intended use.

Do not expose the thermometer to direct sunlight, high temperature/ humidity or any extreme environment, otherwise the function will be reduced.

When the ambient temperature of the thermometer changes too much, such as moving the thermometer from one place of lower temperature to another place of higher temperature, allow the thermometer to remain in a room for 30 minutes where the temperature is between 10°C to 40°C.

8 - DISPOSAL

1) Used batteries should not be disposed of in the household rubbish.

Used Batteries should be deposited at a collection point.

2) At the end of its life, the appliance should not be disposed of in household rubbish. Enquire about the options for environment-friendly and appropriate disposal. Take local regulations into account.

9 - WARRANTY

Our company warrants the Mini Multifunction thermometer at the time of its original purchase and for the subsequence time period of 24 months.

The warranty does not cover the followings:

- The device series number label is torn off or cannot be recognized.
- Damage to the device resulting from misconnection with other devices.
- Damage to the device resulting from accidents.
- Changes performed by users without the prior written authorization of the company.
- Batteries and packaging are not covered under warranty.

10 - CALIBRATION

The thermometer is initially calibrated at the time of manufacture. If this thermometer is used according to the use instructions, periodic re-adjustment is not required. If at any time you question the accuracy of temperature measurements, please contact us timely.

11 - EMC DECLARATION

This equipment needs to be installed and put into service in accordance with the information provided in the user manual.

This product needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided, and this unit can be affected by portable and mobile RF communications equipment.

- Do not use a mobile phone or other devices that emit electromagnetic fields, near the unit. This may result in incorrect operation of the unit.
- Warning: This unit has been thoroughly tested and inspected to assure proper performance and operation!
- Warning : Use of this equipment adjacent to or stacked with other equipments should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

Guidance and manufacturer's declaration - electromagnetic emission

The Mini multifunction thermometer is intended for use in the electromagnetic environment specified below. The customer or the user of the Mini multifunction thermometer should assure that it is used in such an environment.

Emission test	Compliance
RF emissions CISPR 11	Group 1
RF emission CISPR 11	Class B
Harmonic emissions IEC 61000-3-2	Not applicable
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable

Guidance and manufacturer's declaration - electromagnetic immunity

The Mini multifunction thermometer is intended for use in the electromagnetic environment specified below. The customer or the user of Mini multifunction thermometer should assure that it is used in such an environment.

Anti-interference detection	IEC 60601 test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	Contact: ± 8 KV Air: ± 2 , ± 4 , ± 8 , ± 15 KV	Contact: ± 8 KV Air: ± 2 , ± 4 , ± 8 , ± 15 KV
Electrical fast transient/burst IEC 61000-4-4	The input a.c. power ports: ± 2 KV The input d.c. power ports: ± 2 KV Signal input/output ports: ± 1 KV	
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Not applicable

Voltage dips IEC 61000-4-11	0.5 cycles for > 95% (sync angle (degrees):0, 45, 90, 135, 180,225, 270, 315) 1 cycles for >95% UT (sync angle (degrees):0) 25 (50Hz)/30 (60Hz) cycles for 30% U_T (sync angle (degrees):0)	Not applicable
Voltage interruption IEC 61000-4-11	250 (50Hz)/300 (60Hz) cycles for >95% U_T (sync angle (degrees):0)	Not applicable
Power frequency (5 0 H z / 6 0 H z) magnetic field IEC 61000-4-8	30A/m	30A/m

NOTE U_T is the a.c. mains voltage prior to application of the test level.

Guidance and manufacture's declaration - electromagnetic immunity

The Mini multifunction thermometer is intended for use in the electromagnetic environment specified below. The customer or the user of the Mini multifunction thermometer should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80MHz	Not applicable
Radiated RF IEC 61000-4-3	Professional healthcare environment: 3 V/m Home healthcare environment: 10 Vm 80 MHz to 2700 MHz	Professional healthcare environment: 3 V/m Home healthcare environment: 10 Vm 80 MHz to 2700 MHz 

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a- Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy.

To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Non Contact Infrared Body Thermometer is used exceeds the applicable RF compliance level above, the Non Contact Infrared Body Thermometer should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Non Contact Infrared Body Thermometer.

b- Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Guidance and manufacturer's declaration - RF wireless communication equipment immunity

Test frequency (MHz)	Band a) (MHz)	Service a)	Modulation b)	Maximum power (W)	Distance (m)	IMMUNITY
385	380-390	TETRA 400	Pulse modulation ^{b)} 18 Hz	1.8	0.3	27
450	430-470	GMRS 460, FRS 460	FM ^{c)} ± 5 kHz deviation 1kHz sine	2	0.3	28
710	704-787	LTE Band 13, 17	Pulse modulation ^{b)} 217 Hz	0.2	0.3	9
745						
780						
810	800-960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation ^{b)} 18 Hz	2	0.3	28
870						
930						
1720	1700-1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation ^{b)} 217 Hz	2	0.3	28
1845						
1970						
2450	2400-2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation ^{b)} 217 Hz	2	0.3	28
5240	5100-5800	WLAN 802.11 a/n	Pulse modulation ^{b)} 217 Hz	0.2	0.3	9
5500						
5785						

a) For some services, only the uplink frequencies are included.

b) The carrier shall be modulated using a 50 % duty cycle square wave signal.

c) As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

10 - BODY TEMPERATURE

For consumer use only, not for professional use.

- Body temperature varies from person to person and fluctuates during the course of the day. For this reason, it is suggested to know one's normal, healthy forehead temperature to correctly determine the temperature.
- Body temperature runs approximately from 35.5°C to 37.8°C. To determine if one has a fever, compare the temperature detected with the person's normal temperature. A rise over the reference body temperature of 1°C or more is generally indication of fever.
- Different measurement sites (rectal, axillary, oral, forehead, auricular) will give different readings. Therefore it is wrong to compare the measurement taken from different sites.
- Here below are typical temperatures for adults, based on different measurement sites:
 - Rectal: 36.6°C to 38°C
 - Axillary: 34.7°C to 37.3°C

ASTM laboratory accuracy requirements in the display range of 37°C for IR thermometers is $\pm 0.2^{\circ}\text{C}$ ($\pm 0.4^{\circ}\text{F}$), whereas for mercury-in-glass and electronic thermometers, the requirement per ASTM Standards E667-86 and E1112-86 is $\pm 0.1^{\circ}\text{C}$.

Caution: This infrared thermometer meets requirements established in ASTM Standard (E1965-98) Except of clause 5.2.2. It displays subject's temperature over a range of 34.0~43.0°C.

Full responsibility for the conformance of this product to the standard is assumed by (Hetaida Technology Co., Ltd. Add: Room 801, 802, 803, 804, 901, 2# Building Scientific Research Center, Songhu Intelligent Valley, No.6 Minfu Road, Liaobu Town, Dongguan City, Guangdong Province, P.R.China.)

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DOREL FRANCE S.A.S. Z.I. - 9 bd du Poitou, BP 905, 49309 Cholet Cedex, FRANCE
Share Info Consultant Service LLC